TRANSMITTAL LETTER TO THE UNITED STATES DESIGNATED/ELECTED OFFICE (DO/EO/US) CONCERNING A FILING UNDER 35 U.S.C. 371 INTERNATIONAL APPLICATION NO. INTERNATIONAL FILING DATE PRIORITY DATE CLAIMED PCT/US00/09493 April 13 1999 April 11 2000 TITLE OF INVENTION Indicators for monitoring sterilization with Plasma APPLICANT(S) FOR DO/EO/US <u>NGorhanbhai</u> PATEL, Applicant herewith submits to the United States Designated/Elected Office (DO/EO/US) the following items and other information: 1. This is a FIRST submission of items concerning a filing under 35 U.S.C. 371. 2. This is a SECOND or SUBSEQUENT submission of items concerning a filing under 35 U.S.C. 371. This is an express request to begin national examination procedures (35 U.S.C. 371(f)). The submission must include items (5), (6), (9) and (21) indicated below. The US has been elected by the expiration of 19 months from the priority date (Article 31). 5. A copy of the International Application as filed (35 U.S.C. 371(c)(2)) is attached hereto (required only if not communicated by the International Bureau). has been communicated by the International Bureau. is not required, as the application was filed in the United States Receiving Office (RO/US). 6. An English language translation of the International Application as filed (35 U.S.C. 371(c)(2)). is attached hereto. has been previously submitted under 35 U.S.C. 154(d)(4). Amendments to the claims of the International Aplication under PCT Article 19 (35 U.S.C. 371(c)(3)) are attached hereto (required only if not communicated by the International Bureau). have been communicated by the International Bureau. have not been made; however, the time limit for making such amendments has NOT expired. have not been made and will not be made. 8. An English language translation of the amendments to the claims under PCT Article 19 (35 U.S.C. 371 (c)(3)). 9. X An oath or declaration of the inventor(s) (35 U.S.C. 371(c)(4)). An English lanugage translation of the annexes of the International Preliminary Examination Report under PCT Article 36 (35 U.S.C. 371(c)(5)). Items 11 to 20 below concern document(s) or information included: 11. An Information Disclosure Statement under 37 CFR 1.97 and 1.98. 12. An assignment document for recording. A separate cover sheet in compliance with 37 CFR 3.28 and 3.31 is included. 13. A FIRST preliminary amendment. 14. A SECOND or SUBSEQUENT preliminary amendment. 15. A substitute specification.

A computer-readable form of the sequence listing in accordance with PCT Rule 13ter.2 and 35 U.S.C. 1.821 - 1.825.

A second copy of the English language translation of the international application under 35 U.S.C. 154(d)(4).

A second copy of the published international application under 35 U.S.C. 154(d)(4).

16.

17.

18.

19.

20. 🗍

A change of power of attorney and/or address letter.

Other items or information:

U.S. APPLICATION NOT GENERAL	4040 °	NTERNATIONAL APPLICATION NO		ATTORNEY'S	DOCKET NUMBER
21. The following	ng fees are submitted:	· · · · · · · · · · · · · · · · · · ·		CALCULATIONS	FTO USE ONLY
Neither international nor international sea	FEE (37 CFR 1.492 (a) I preliminary examination for the (37 CFR 1.445) arch Report not prepare	on fee (37 CFR 1.482)	· · · · · · · · \$1040.00		
International prelimi	nary examination fee (3	37 CFR 1.482) not paid to epared by the EPO or JPO)		
International preliming but international sear	nary examination fee (3 rch fee (37 CFR 1.445(a	37 CFR 1.482) not paid to a)(2)) paid to USPTO	USPTO \$740.00		
International preliming but all claims did not	nary examination fee (3) satisfy provisions of PC	37 CFR 1.482) paid to US CT Article 33(1)-(4)	SPTO \$710.00		
and all claims satisfie	ed provisions of PCT A	37 CFR 1.482) paid to US rticle 33(1)-(4)	\$100.00	\$ 890	
Surcharge of \$130.00 nonths from the earlie	for furnishing the oath est claimed priority date	or declaration later than e (37 CFR 1.492(e)).	<u>20</u> 30	\$ 130	
CLAIMS	NUMBER FILED	NUMBER EXTRA	RATE	\$	
Total claims	41 - 20 =	21	x \$18.00	\$ 378	
ndependent claims	- 3 =		x \$84.00	\$	
ULTIPLE DEPEND	ENT CLAIM(S) (if app	·	+ \$280.00	\$	
		OF ABOVE CALCU		\$ 1398	
Applicant claims are reduced by 1/		e 37 CFR 1.27. The fees	indicated above +	\$	}
		<u> </u>	JBTOTAL =	699	
rocessing fee of \$130 conths from the earlie	0.00 for furnishing the I est claimed priority date	English translation later th		\$	
		TOTAL NATIO	NAL FEE =	\$ 699	
ee for recording the eccompanied by an app	enclosed assignment (37) propriate cover sheet (3	7 CFR 1.21(h)). The assig 7 CFR 3.28, 3.31). \$40.0	gnment must be 00 per property +	\$	
		TOTAL FEES EN	NCLOSED =	\$ 699	
				Amount to be refunded:	\$
				charged:	\$
A check in th	ne amount of \$	to cover the	e ahove fees is enclos	e d	
o. 🔯 Please charge	e my Deposit Account I copy of this sheet is enci	No. 19-1218 in	the amount of \$	to cover	the above fees.
			ate copy of this sheet	is cholosed.	
overpayment d. Fees are to be	e charged to a credit ca	rd. WARNING: Informa on this form. Provide cr	ation on this form ma	y become public. (Credit card on PTO-2038.
overpayment d. Fees are to be information s NOTE: Where an a	e charged to a credit ca should not be included ppropriate time limit i	rd. WARNING: Informa	ation on this form may redit card information	y become public. (and authorization	on PTO-2038.
overpayment Fees are to be information s NOTE: Where an ap 1.137 (a) or (b)) mus	e charged to a credit ca should not be included ppropriate time limit to the filed and granted	rd. WARNING: Informa on this form. Provide cr	ation on this form may redit card information	y become public. (and authorization	on PTO-2038.
overpayment Fees are to be information s NOTE: Where an ap 1.137 (a) or (b)) mus	e charged to a credit ca should not be included ppropriate time limit u it be filed and granted IDENCE TO:	rd. WARNING: Informa on this form. Provide cr	ation on this form may redit card information	y become public. (and authorization et, a petition to re	on PTO-2038.
overpayment Fees are to be information s NOTE: Where an application of the control of the cont	e charged to a credit can chould not be included ppropriate time limit until be filed and granted EDENCE TO: hr kince Kim PC	rd. WARNING: Informa on this form. Provide cr	ation on this form may redit card information .495 has not been m on to pending status.	y become public. (and authorization et, a petition to re	on PTO-2038.
overpayment Fees are to be information s NOTE: Where an all 1.137 (a) or (b)) mus SEND ALL CORRESPON Omri M. Bell Calitto, Bell PO Box 1477	e charged to a credit can chould not be included ppropriate time limit us to be filed and granted EDENCE TO: hr & Kim PC	rd. WARNING: Information this form. Provide or under 37 CFR 1.494 or 1 to restore the application	ation on this form mare dit card information. 495 has not been much to pending status.	y become public. (and authorization et, a petition to re	on PTO-2038.
overpayment d. Fees are to be information s NOTE: Where an application of the control of the c	e charged to a credit can chould not be included ppropriate time limit us to be filed and granted EDENCE TO: hr & Kim PC	rd. WARNING: Informa on this form. Provide cr	ation on this form mare dit card information. 495 has not been much to pending status. SIGNATUR NAME ON REGISTRA	y become public. (and authorization et, a petition to re	on PTO-2038.
overpayment Fees are to be information s NOTE: Where an ap 1.137 (a) or (b)) mus SEND ALL CORRESPON Omri M. Bel Octivto, Bel PO Box 1477 Edison NJ (e charged to a credit ca should not be included ppropriate time limit us at be filed and granted NDENCE TO: hr hr & Kim PC 7	rd. WARNING: Information this form. Provide crumder 37 CFR 1.494 or 1 to restore the application of the state	ation on this form mare dit card information. 495 has not been much to pending status. SIGNATUR NAME ON REGISTRA	y become public. (and authorization et, a petition to re	on PTO-2038.
overpayment i. Fees are to be information s NOTE: Where an application of the control of the c	e charged to a credit ca should not be included ppropriate time limit us at be filed and granted NDENCE TO: hr hr & Kim PC 7	rd. WARNING: Information this form. Provide crumder 37 CFR 1.494 or 1 to restore the application of the state	ation on this form mare dit card information. 495 has not been much to pending status. SIGNATUR NAME ON REGISTRA	y become public. (and authorization et, a petition to re	on PTO-2038.
overpayment Fees are to be information s NOTE: Where an application of the control of the cont	e charged to a credit ca should not be included ppropriate time limit us at be filed and granted NDENCE TO: hr hr & Kim PC 7	rd. WARNING: Information this form. Provide crumder 37 CFR 1.494 or 1 to restore the application of the state	ation on this form mare dit card information. 495 has not been much to pending status. SIGNATUR NAME ON REGISTRA	y become public. (and authorization et, a petition to re	on PTO-2038.
overpayment Fees are to be information s NOTE: Where an application or (b)) must send all correspondent of the control of th	e charged to a credit ca should not be included ppropriate time limit us at be filed and granted NDENCE TO: hr hr & Kim PC 7	rd. WARNING: Information this form. Provide crumder 37 CFR 1.494 or 1 to restore the application of the state	ation on this form mare dit card information. 495 has not been much to pending status. SIGNATUR NAME ON REGISTRA	y become public. (and authorization et, a petition to re	on PTO-2038.
overpayment Fees are to be information s NOTE: Where an ap. 137 (a) or (b)) mus END ALL CORRESPON Omri M. Bel PO Box 1477 Edison NJ (e charged to a credit ca should not be included ppropriate time limit us at be filed and granted NDENCE TO: hr hr & Kim PC 7	rd. WARNING: Information this form. Provide crumder 37 CFR 1.494 or 1 to restore the application of the state	ation on this form mare dit card information. 495 has not been much to pending status. SIGNATUR NAME ON REGISTRA	y become public. (and authorization et, a petition to re	on PTO-2038.

TITLE OF THE INVENTION

Indicators for monitoring sterilization with plasma

10

BACKGROUND OF THE INVENTION 1. FIELD OF THE INVENTION

The present invention relates to chemical indicators for monitoring oxidizing vapors or plasmas, in particular for sterilization of medical supplies with plasma and detection of hydrogen peroxide.

15

2. BRIEF DESCRIPTION OF PRIOR ART

A wide variety of medical supplies and biological wastes are sterilized with materials and techniques, such as steam, ethylene oxide (ETO), high energy radiation and plasma. It is essential to assure that these supplies or wastes are actually sterilized. A number of indicators, dosimeters and monitors have been proposed in the literature. They include biological and chemical indicators. The color changing chemical indicators are inexpensive and hence are widely used.

25

30

35

20

In order to assure the sterilization with plasma has taken place, the indicator must determine integral value of parameters, such as time, temperature, humidity, and concentration of the plasma gas. Biological indicators made from cultures, such as Bacillus subtilis spores, bacillus pumilus spores and clostridium sporogenes spores have been used for monitoring the sterilization, for example, US patents 5,801,010, 5,788,925, and 5,866,356. However, chemical indicators are preferred because they are simple and inexpensive.

W09846994A1 and W09846279A1 describe compositions comprising dyes and organic amines which change color when contacted with hydrogen

5 peroxide plasma. W09852621A1 describes dyes, such as, acid fuschin which change color when contacted with plasma of hydrogen peroxide.

SUMMARY OF THE INVENTION

There is provided a polymer device for monitoring plasmas as well as the presence of oxidants, comprising:

at least one layer of polymer, having incorporated therein

- a) an indicator capable of undergoing at least one color change
- b) an activator for said indicator wherein said activator, when contacted with said plasma or oxidant, introduces a reaction wherein the product of said reaction causes said indicator to undergo said color change.

The indicators suitable for use in this device include pigments, dyes, precursors of said dyes, and mixtures of any of these. Suitable indicators include pH-sensitive sensitive dyes such as phenol red, m-cresol purple, pararosaniline or mixtures thereof. A desirable quality of the indicator is the ability to undergo halogenation, especially bromination or oxidation. The source of such brominating or oxidizing agents being the reaction between the plasma and the activator. Desirably the indicator undergoes a yellow-to-blue, red-to-yellow or red-to-blue color change.

25

30

35

15

20

The polymer used in the device is, suitably, soluble in water or dispersible in an aqueous medium solvent. A broad class of polymers may be used. They may be homopolymers, copolymers or a mixture thereof, suitably a polymer of styrene, acrylate, acrylic acid, acrylamide, vinyl acetate, vinyl alcohol, vinyl chloride, styrene, polyurethanes, cellulose nitrate, carboxymethyl cellulose or a mixture thereof. Desirably, the polymer is an acrylate polymer, cellulose nitrate or carboxymethylcellulose

Suitably the reaction product of the activator and the plasma is a haloacid. Suitable activators are salts such as halides, preferably bromides. Desirable halides include bromides of alkali metals or of quaternary amines such as tetrabutylammonium bromide or tetraethylammonium bromide or a mixture thereof. The activator may also be a salt of an amine and an organic or inorganic acid. It may also be a metal thiocyanate such as sodium thiocyanate.

10 . The device may additionally comprise an additive to control the diffusion of plasma gases, such a crosslinking agent or a plasticizer. These may include zinc compounds or polyaziridines

The device may have two layers, that is to say additionally comprising a polymeric top layer. This may be a wedge shaped polymeric top layer.

The process of making a device of the present invention comprises dissolving the components thereof in a solvent thereof, applying the thus formed solution to a substrate and permitting the solvent to evaporate.

20

The substrate may be a container for an item to be sterilized. It may also be a plastic film, paper or metal, including but not limited polyester film or spun bonded polyolefins.

In a desirable embodiment of the invention the solution is an ink formulation suitably an aqueous ink formulation most suitably one which comprises an acrylate polymer.

A process of using a device of the present invention for monitoring sterilization of materials comprises the steps of affixing the device to said materials or containers containing same, carrying out the process of sterilization including the step of introducing the plasma into a vessel containing said materials or containers thereof and observing the presence of a color change of said device.

Among the plasmas which may be utilized in this process are those derived from hydrogen peroxide, perchloric acid and oxygen, most suitably from hydrogen peroxide.

DETAILED DESCRIPTION OF THE DRAWINGS

Figure 1. A side schematic cross section of one embodiment of the plasma sterilization indicator of the invention where an indicator layer comprised of a polymeric binder, plasma activator and plasma indicator is applied on a substrate.

15

10

Figure 2. A side schematic cross section of the plasma sterilization indicator of the invention having an adhesive layer and a release layer.

Figure 3. A side cross section of a two layered plasma sterilization indicator 20 layer.

Figure 4. A side schematic cross section of the moving boundary plasma sterilization indicator created by coating a wedge shaped barrier on the plasma indicator.

25

30

35

These devices can also be used for detecting and monitoring oxidants, suitably in gaseous form, such as hydrogen peroxide ozone, and the like.

DESCRIPTION OF THE PREFERRED EMBODIMENTS

The device can be best described by reference to the Figures. As shown in Figure 1, the device in its simplest form is comprised of an indicator layer 20, applied on a substrate 10. The substrate 10 can also be a container, such as pouch for products to be sterilized. The indicator layer 20 is composed of a polymer 50 containing at least one indicator 30 capable of undergoing a color change when reacted with a reactive species produced by

exposure of the activator <u>40</u> to plasma. The layer <u>20</u> may additionally contain other additives <u>80</u> to control the rate of the color change, such as a crosslinking agents or plasticizers to control diffusion of plasma gases either to increase or decrease the penetration of plasma.

As shown in Figure 2, the substrate <u>10</u> of the device can further comprise an adhesive layer <u>60</u>. The adhesive layer allows the device to be affixed to a container of product to be sterilized. To the bottom of the adhesive layer <u>60</u>, can further be affixed a release layer <u>70</u> for ease in packaging and for removal just prior to use.

15

20

10

As shown in Figure 3, a simple form of a multi-layer device can have two indicator layers <u>20</u> and <u>200</u>. The multi-layer device can thus undergo more than one color change. The device can be created by applying another indicator layer <u>200</u> over the first indicator layer <u>20</u>. The indicator <u>300</u> of the indicator layer <u>200</u> would be different from indicator <u>30</u> of the indicator layer <u>20</u>. The polymers <u>50</u> and <u>500</u>, additives <u>80</u> and <u>800</u> and the activators <u>40</u> and <u>400</u> of the first indicator layer <u>20</u> and the second indicator layer <u>200</u> could be the same or different. When exposed to plasma, the top indicator layer <u>200</u> may undergo a color change (e.g., red-to-colorless) first followed by the color change of the bottom indicator layer <u>20</u> (e.g., colorless-to-blue). Thus the device would undergo two color changes. There are many variations of the multi-layer device. For example, the device can have more than two indicator layers and can have one or more barrier layers between or on them.

30

35

25

A simple form of a moving boundary indicator device is shown in Figure 4. The moving boundary device can be created by applying a wedge shaped barrier layer 2000 over the indicator layer 20. The barrier layer will resist but will be permeable to plasma and its components. There are many variations to this moving boundary device. For example, the wedge barrier coat may contain indicator and both the indicator and the barrier layers be in form of wedge but in opposite direction.

10

15

Other variations of the plasma indicator device are also possible, for example, a gradient device can be created by coating a series of formulations having the time required for the color change either increase or decrease. Such gradient can be obtained by coating such formulations in form of lines or bars next to each other.

The device may also include a printed number, an image, a bar code or a message, e.g., "if this print is green, the product inside is sterilized".

Activators produce other reactive species, such as those from compounds having: -I, -Br, -Cl, -NH-, -NH₂, -OH, =O, -NO, -NO₂, -NO₃, -SO₃, -SO₃H, -COOH, -PO₄, SCN, etc. Especially desirable are bromides, i.e., from tetraethyl ammonium bromide (TEAB), add Br as an oxochrome to the dyes. Similarly other activators add other oxochromes to the dye.

20

25

30

35

Results from a series of experiments indicate the product formed by the reaction of phenol red and the bromide salts is most likely bromophenol blue. We have also found that nature, e.g., pH, of the medium i.e., binder, plays an important role in determination of color of the reaction product, especially the dye produced is a pH-sensitive dye. For example, when exposed to hydrogen peroxide, phenol red changes from yellow-to-orange red in polyvinyl alcohol while it changes from yellow-to-blue in EC001270 (an acrylate printing ink extender supplied by Environmental Inks and Coating Co., Lithicum, MD). The pH of EC001270 coating is about 4-5 while that of polyvinyl alcohol is about 7. Bromophenol blue appears red in polyvinyl alcohol while blue in EC001270.

The color change in EC001270 and TEAB system for chlorophenol red is yellow-to-purple, m-cresol purple is yellow-to-blue, cresol red is yellow-to-purple, phenol red is yellow-to-blue, and thymol blue is yellow-to-blue upon exposure to hydrogen peroxide vapor and its plasma. Evaluation of these and

7

other dyes which develop colors indicates that brominated dyes are formed for some of these dyes upon exposure to hydrogen peroxide and its plasma.

As shown in Table 1 below, the conclusion that EC001270, TEAB and dyes which develop color is supported by the pH range and colors of the substrate dyes (dyes before the treatment) and corresponding brominated dyes (after exposing to hydrogen peroxide and its plasma).

Table 1. Color of the starting substrate dyes and corresponding brominated dyes at low and high pH ranges.

Substrate dye	Low	Color	High PH	Color	High Color Brominated dye pH	Low color pH		High PH	Color
for locadacard	48	Yellow 6.4	6.4	Red	Bromochlorophenol red	3.2	yellow 4.8	4.8	Purple
		Yellow	0.6			3.8	yellow	5.4	Blue-green
Creed red	7.2	Yellow		Purpl	Purple Bromocresol purple	5.2	yellow	8 .	6.8 Purple
Cleson ca Phenol red	8	Yellow		Red	Bromophenol blue	3.0	yellow	4.6	Bine
Thymol blue	8.0	Yellow	9.2	Blue	Bromothymol blue	6.0	yellow	9.7	Blue
Xvlenol blue	8.0	Yellow	9.6	Blue	Bromoxylenol blue	6.0	yellow		7.6 Blue

<u>Indicators:</u> Any material which undergoes a color change in presence of a plasma activator can be used a plasma indicator. Plasma indicators are also referred herein to as indicators. Most preferred classes of plasma indicators are dyes, pigments and their precursors. A large number of dyes and dye precursors, e.g., are listed in Table 2.

10

5

Table 2. List of useable dyes.

These have been tested with tetrabutylammonium bromide as an activator in polyacrylate (EC001270) and in the presence of solvent based cellulose nitrate.

15

20

25

30

35

Acid alizarin violet N, acid black 24, acid black 48, acid blue 113, acid blue 120, acid blue 129, acid blue 161, acid blue 25, acid blue 29, acid blue 40, acid blue 41, acid blue 45, acid blue 80, acid blue 93, acid fuschin, acid green 25, acid green 27, acid green 41, acid orange 74, acid red 1, acid red 114, acid red 151, acid red 88, acid violet 17, acid violet 7, acid yellow 99, acridine orange, acridine orange base, acridine orange G, acridine yellow G, acriflavine hydrochloride, alcian blue 8GX, alcian yellow, alizarin, alizarin blue black SN, alizarin complexone, alizarin complexone dihydrate, alizarin red, alizarin violet 3R, alizarin yellow GG, alizarin yellow R, alkali blue 6B, alkali fast green 10GA, alphazurine A, aluminon, aminoacridine hydrochloride, aminoanthraquinone, aminophthalhydrazide, aniline blue, astra blue 6GLL, auramine O, azocarmine, azocarmine B, azure A, azure B, azure B thiocyanate, azure C, basic blue 3, basic blue 41, basic blue 66, basic fuchsin, basic red 29, basic yellow 11, benzo purpurin 4B, biebrich scarlet NA salt, bismarck brown B, bismarck brown Y, blue tetrazolium, bordeaux R, brilliant blue B, brilliant blue G, brilliant cresyl blue ALD, brilliant crocein MOO, brilliant green, brilliant sulphaflavine, brilliant yellow, bromochlorophenol blue, bromocresol green, bromocresol purple, bromophenol blue, bromopyrogallol red, bromothymol blue, bromoxylenol blue, calmagite, carbol fuchsin, carminic acid, carotene, celestine blue, Chicago sky blue, chlorophenol red, chrome azurol S, chromotrope 2B, chromotrope 2R, chromoxane cyanine B,

10

15

20

25

30

35

chrysoidin, chrysophenine, cibacron brilliant red 3BA, Congo red, copper(II) cresolphthalein, phthalocyanine, cresol purple, cresol red, cresol, cresolphthalein complexone, crystal violet, red. curcumin, darrow RC, dibromofluorescein, diaminoacridine hemisulfate, diazo red dichloroindophenol, dicinnamalactone, dichlorofluorescein, diethyloxacarbocyanine iodide, diethylaminomethyl coumarin. diethylthiatricarbocyanine iodide, dihydroxy benzenesulfonic acid, dilithium blue. dimethylalyoxime. phthalocyanine, dimethyl methylene dimethylindoaniline, dinitro diphenylamine, diphenylthiocarbazone, direct blue 71, direct green 6, direct red 23, direct red 75, direct red 81, direct violet 51, direct yellow 62, disodium phthalocyanine, disperse blue 14, disperse blue 14, disperse blue 3, disperse orange, disperse orange 11, disperse orange 25, disperse yellow 7, emodin, eosin B, eosin Y, eriochrome black T, eriochrome blue black B, erioglaucine, erythrosin B, ethyl eosin, ethyl orange, ethyl red, ethyl violet, Evans blue, fast black, fast blue B salt, fast blue BB, fast blue RR, fast blue RR salt, fast corinth V salt, fast garnet GBC base, fast green FCF, fast red aluminum salt, fast red violet LB salt, fast violet B salt, fat brown RR fat green GDC salt, flavazin I, fluorescein, fluorexon, gallocyanine, guinea green B, hematoxylin, hydroxy naphthol blue, 1.4-hydroxynaphthoguinone, indigo, indigo carmine, indoline blue, iron(II) phthalocyanine, janus green B. lacmoid, leishman stain, leuco crystal violet, leucomalachite green, leucoquinizarin, light green SF yellowish, lissamine green B, litmus, base, malachite green hydrochloride, luxol fast blue, malachite green malachite green oxalate, metanill yellow, methyl eosin, methyl green, methyl orange, methyl red, methyl violet 2B, methyl violet B base, methyl yellow, methylene blue, methylene green, methylene violet 3RAX, methylesculetin, methylthymol blue, mordant blue 9, mordant brown 24, mordant brown 4, mordant orange, mordant orange 1, mordant orange 6, mordant red 19, mordant yellow 10, morin hydrate, murexide, naphthochrome green, naphthol AS, naphthol blue black, naphthol green B, naphthol yellow, naphtholbenzein, naphtholbenzene, naphtholphthalein, neutral red, new coccine, new fuchsin, new methylene blue N, nigrosin, nile blue A, Nile blue

10

15

20

25

30

35

chloride, nitrazine yellow, nitro red, nitro-phenanthroline, nitrophenol-2, nitrophenol-3, nitrophenol-4, nitrophenylazo-resorcinol, nuclear fast red, oil blue N, oil red EGN, oil red O, orange G, orange II, palatine chrome black 6BN, palatine fast yellow BLN, pararosaniline acetate, pararosaniline base, pararosaniline chloride, patent blue VF, pentamethoxytriphenylmethanol, phenanthroline, phenazine, phenol red, phenolphthalein, phenolphthalein diphosphate, phenothiazine, phenylazoaniline, phenylazodiphenylamine, phloxine phenylazoformic acid. phenylazophenol. B. phthalocynine, pinacyanol chloride, plasmocorinth, ponceau S, primuline, procion red MX-5B, procion yellow H-E3G, prussian blue, purpurin, pyridlazo naphthol, pyridylazoresorcinol sodium salt, pyrocatechol violet, pyrogallol red, pyronin B, quinaldine red, quinizarin, quinoline yellow, reactive black 5, reactive blue 15, reactive blue 2, reactive blue 4, reactive orange 16, resazurin, resorcin crystal violet, rhodamine B, rhodamine B base, rhodamine GG, rhodamine S, rhodanine, rosalic acid, rose bengal rose bengal iactone, safranine O, solvent blue 35, solvent blue 59, solvent green 3, styryl 7, sudan black B, sudan orange G, sudan red 7B, sulfobromophthalein sodium salt, tartrazine. tetrabromophenol blue. tetrabromo sulforhodamine B. phenolphthalein, tetrabromo phenolphthalein, tetraiodo phenolphthalein, tetraphenyl-butadiene, tetrazolium violet, thiazol yellow G, thioflavin S, thioflavin T, thionin, thymol blue, thymolphthalein, thymolphthalein monophosphate, thymolphthalein monophosphate, toluidine blue triphenylmethyl bromide, tropaelin O, trypan blue, turmeric, vanillin azine, variamine blue RT salt, variamine blue RT salt, victoria blue B, victoria blue B, victoria pure blue BO, wright stain, xilidine ponceau 2R,, xylenol blue, and xylenol orange.

Some dyes, e.g., malachite green, reduced with ascorbic acid, sodium sulfite and formaldehyde were also used. Some of these dyes are fluorescence dyes and they either lost their fluorescence or there was a change in fluorescence.

10

15

20

and a first party party party from the first party

<u>Simultaneous reactions</u>: We have also observed that certain dyes, e.g., phenol red show development of color (due to formation of bromophenol blue) followed by de-coloration upon prolonged exposure to hydrogen peroxide or plasma treatment.

Activators: Any chemical which produces reactive species with a plasma or hydrogen peroxide which, when reacted with a indicator compound, degrade, cleave or add to the molecule or react to form a colored compound can be used as an indicator activator. Similarly, a chemical which when reacted with plasma or hydrogen peroxide and produces a chemical which changes environment, e.g., pH of the medium, can also be used as an indicator activator. Production of an acid or base introduces a change of a pH-sensitive dye. A variety of classes of organic and inorganic compounds may be used as activators for monitoring hydrogen peroxide and its plasma. They include alcohols, amides, amines, bisulfites, bisulfates, carbonates, carbamates, chelates, metal complexes, cyanates, esters, halocarbons, ketones, nitrites, nitrates, nitriles, nitro, nitroso, oximes, phenols, phosphates, sulfates, sulfides, sulfites, thiocyanates, ureas, urethanes, salts, oxidants and reducing agents. Organic and inorganic salts, especially halides were very effective activators. The specific examples of compounds explored as plasma activators with some selected dyes (e.g., chrome azurole S, methyl green, acid fuschin, direct blue 71 and nile blue) are listed in Table 3.

Table 3. List of specific examples of activators.

30

35

25

Abietic acid, acetone oxime, aluminum acetylacetonate, aluminum ammonium sulfate, aluminum chloride, aluminum sulfate, amino deoxy d-sorbitol, ammonium acetate, ammonium bisulfite, ammonium bromide, ammonium carbamate, ammonium nitrate, ammonium sulfamate, ammonium sulfite, ammonium thiocyanate, ammonium thiosulfate, ascorbic acid, azodicarbonamide, azodicarbonamide, benzilic acid, benzoic acid,

5 acid, benzophenone, benzophenone tetracarboxylic benzophenonetetracarboxylic benzoguinone dioxime. diahydride, benzoquinone dioxime, benzyloxy)phenol, butyl phenol, caffeine, calcium ferrocyanide, catechol, catechol, chloranilic acid, copper thiocyanate, cupferron, cupferron, cyclopenatanone oxime, dehydroacetic acid, di-butyl-t-10 4-methylphenol, dihydroxy acetophenone, dihydroxy dimethoxy disulfonic acid, benzophenone, dihydroxy naphthalein dihydroxyacetophenone, dihydroxy-dimethoxybenzophenone, dimethyl fumarate, dimethyl tartrate, diphenyl butyro lactone, diphenylglyoxime, dithizone or diphenylthiocarbazone, di-t-butyl-4-methylphenol, diphenylthiocarbazone, ethylcarbonate, ethylenediamine tetraaceticacid and 15 ts salts, ferroin, fumaric acid, gallic acid, gluconic acid fe(ii) salt, glucose penta acetate, glucose penta acetate, glucose pentaacetate, glutaaric acid, glycerophosphate, glyconolactone, hexahloro norborene dicarboxylic acid, hydroquinone, hydroxy acetophenone, hydroxy acetophenone, hydroxy hydroxy methoxybenzophenone, octvloxy 20 cinnamic acid, hydroxy hydroxybenzophenone, benzophenone, hydroxybenzophenone, hydroxymethoxybenzophenone, hydroxyquinoline, hydroxyquinoline, inositol, iron acetylacetonate, iron complexes such as potassium ferrocyanide, iron sulfate, isoascorbic acid, levulinic acid, maleic acid, maleic acid, malic acid, 25 mandelic mercaptobenzothiazole. methyldinitrosalicilate, acid. methyldinitrosalicylate, methylesculetin, methyltrihydroxybenzoate, naphthol, naphthol-disulfonic acid, naphthoquinone tetrasulfate sodium salt, nitron, nitron, nitroso-1,2-naphthol, nitrosophenol, oxalic acid, phenanthroline, phenanthroline, phthalide, propylgalliate, propylgalliate, pydine aldoxime, pyruvic acid, resorcinol, rutin hydrate, salicyladoxime, salicyladoxime, 30 salicylanamid, salicylanilide, salicylic acid, sodium acetylacetonate, sodium sodium cyante. sodium diethyldithio carbamate, bisulfite. diethyldithiocarbamate, sodium dithionite, sodium hydrosulfide, sodium nitrite, sodium persulfate, sodium sulfite, sodium sulfite, sodium sulfite, sodium 35 sodium thiocyanate, sodium thiocyanate, sodium thiosulfate, sulfite. sulfosalicyclic acid 5, tannic acid, tetrabutylphosphonium bromide,

tetrahydroxybenzophenone, tetrahydroxybenzophenone, tetramethylhexane diamine, tetronic acid, tetronic acid, thiodiglycolic acid, thiodipropionic acid, thioglycolic acid, thiourea, tribenzylamine, trichloroacetamide, trichlorobenzylacetate, trihydroxybenzophenone, trihydroxybenzophenone, urea, vitamin-c, and vitamin-c palmitate

10

Effect of halides: Many organic and inorganic halides are highly sensitive indicator activators. These halides include, acetyl choline chloride, ammonium bromide, choline chloride, choline iodide, dodecyltrimethylammonium bromide, glycidil trimethyl ammonium chloride, potassium bromide, potassium iodide, sodium iodide, tetrabutyl ammonium iodide, tetraethyl ammonium bromide, tetrahexyl ammonium bromide, tetramethyl ammonium chloride and tetrabutyl phosphonium bromide. Bromides were more effective than other halides.

20

15

Acid base salts, e.g., those produced by neutralizing, primary, secondary and tertiary amines (e.g., hexyl amine, diethanolamine, tetramethylhexane diamine) with acids such as acetic acid, hydrochloric acid, hydrobromic acid and hydroiodic acid are also effective activators.

25

The time required for the color change can be controlled by using a proper mixture of the activators/halides/bromides. In order to minimize effect of ethylene oxide lower concentration of halides is preferred for certain dyes.

30

Effect of concentration of activator: With certain dyes as low as 0.5% of activator is effective in introducing noticeable color change with plasma. The rate of color change increases with increasing the concentration of the activator. High concentrations of activator sometimes lead to decoloration of the colored formed. 0.5 to 50 w/w% concentration of a activator may be useful. The preferred concentration range is 2-10 w/w% of the total solid.

10

15

20

25

30

35

Polymers: A matrix or a medium in which a wide variety of organic and inorganic polymeric materials can be used as carrier for the indicator as long as the activators and indicators can be dissolved or dispersed in them. Water-soluble, water-dispersible and water-insoluble polymers may be used. It is preferable to use water-soluble and water-dispersible polymers. These can be formulated in the form of ink formulations, such as flexo and gravure inks. Other inks such as those for letter press, offset and screen printing can also be used. Selection of a polymer depends upon the printing/coating equipment to be used.

Usually acrylic polymers, emulsion of acrylic polymers, occasionally natural polymers, such as starches, lignins, and lignin derivatives are used in inks. Resins are water soluble or emulsifiable through neutralization with basic compounds, such as ammonia and amines. Inks contain a variety of additives to eliminate foaming, dispersion of pigments, rheological modifiers, and slip agents.

Polymers used in inks include homopolymers, copolymers and block-copolymers including those of ethylene acrylic acid, ethylene methacrylic acid, ethylene n-butyl acrylate, and ethylene methyl acrylate. The polymers used in the inks could also be a mixture of homo and copolymers, e.g., those of methylmethacrylate, acrylic acid, styrene, methyl acrylate, other esters may include crosslinking agents, such as polyvinylazaridines metal salts and complexes, e.g., that of zinc.

Commercial sources for suitable polymers for ink formulations include Air products (Allentown, PA), Rohm and Haas (Philadelphia, PA), S.C. Johnsons and Sons (Racine, WI), Witco (Houston, PA) and ESI (Valley Stream, NY). Though a large number of polymers are suitable an ink extender, EC001270 made by Environmental Inks and Coating Co., Lithicam, MD which is composed about 40% styrene-acrylic polymers, a few percent

15

20

25

30

35

5 ammonium hydroxide, additives, such as wax and alcohol and the balance water, has been found very suitable.

The nature of the polymer, e.g., pH and permeability to the gases, plays an important role. Color change of phenol red and TEAB system with hydrogen peroxide vapor is faster in EC001270 (an acrylate) and slower in Witcobond 213, (a polyurethane of Witco Corporation, Houston, TX). Permeability of the reactive gases can also be controlled by addition of a crosslinking agent and other additives such as plasticizers which change the rate of diffusion. When exposed to hydrogen peroxide vapor, the color change of phenol red with potassium iodide in EC001270 is blue while it is red in polyvinyl alcohol. The pH of the coating of EC001270 is about 4-5 while that of polyvinyl alcohol is about 7.

Though aqueous ink or coating formulations are preferred, one can use solvent based coating formulations polymers used in such formulations are cellulose nitrate, carboxymethyl cellulose, polyolefins, polyvinylchloride, polyurethane, polysilicones and polyepoxy.

Effect of nature of the polymer: Coatings prepared from phenol blue as indicator, TEAB as activator and a variety of polymeric materials as binders were exposed to hydrogen peroxide vapor. EC001270 extender was the fastest binder, Witco 160 (a polyurethane of Witco Corporation) and Rovace 571 (a polyvinylacetate latex by Rohm and Haas) were medium and Airflex TL40 (a polyvinylacetate latex by Air Products) was the slowest. They showed the intermediate color, purple before turning blue.

If the barrier coat is of varied thickness, a moving boundary device can be created. Such device was created by applying a wedge shaped coating with a wedge shaped coating bar of EC001270 on the coating of EC001270, phenol red and TEAB. A boundary was created between the original yellow color and the new blue color after about 2 hours when exposed to hydrogen

5 peroxide. The boundary moved from the thin end towards the thick end with prolonged exposure.

Two-layer device: The device may have more than one indicator layer, each containing indicator, activator and polymer. In order to get at least more than one color change the indicator should be different in different indicator layers and should undergo different color changes.

Both layers do not have to undergo color changes with plasma. Even if one layer undergoes a color change, a color change can be noticed, especially if the top layer becomes colorless

<u>Mixtures</u>: Desired colors, color changes and time required for the color changes can be obtained by mixing proper indicators, activators, additives and polymer in appropriate amounts.

20

25

30

35

15

Multi color changes: We have found that by mixing more than one dye, each undergoing different color changes, one can obtain more than one color changes. Similarly, one can use a mixture of pigments or dyes which are not affected by plasma and indicators, especially those which become colorless or lighter color with plasma, to get more than one color change. The rate of reaction was also varied by using a mixture of more than one activators of different reactivities.

<u>Substrate</u>: Though the device may be self-supporting polymer film containing the activator and indicator, it is desirable to prepare the device on substrate. The device can be made by coating the indicating formulation on a substrate. The substrate could be any solid surface; for example one made from paper, plastic, ceramic or metal. The substrate could be a container, e.g., bag or pouch, for items to be sterilized. The sterilization indicator can also be prepared in form of stickers, tapes and strips.

10

15

20

25

30

35

Although any solid substrate having a smooth surface can be used, a preferred substrate is a flexible and transparent plastic film, natural (cellulose) or synthetic (e.g., spun bonded polyolefins, e.g., Tyvek^R) papers. Plastic films, such as polyethylene, polypropylene, polyvinylchloride, polymethylmethacrylate, polyurethanes, nylons, polyesters, polycarbonates, polyvinylacetate, cellophane and esters of cellulose can be used as the transparent substrate. Metal foils, such as aluminum can be used. The most preferred substrates are the 5 - 300 microns thick films of polyethylene terephthalate and Tyvek^R.

<u>UV absorber</u>: If a plasma indicator formulation is sensitive to ultraviolet light, UV absorbing materials can be added to minimize the effect. A large number of UV absorbers are available commercially, e.g., those used in sun-tan lotions, e.g., Tinuvins of Ciba-Geigy Corporation. UV absorbers include compounds such as maleic acid, sodium salicylate, benzophenone, or benzophenone tetracarboxylate and a large number of polyaromatic compounds.

<u>Temperature</u>: The indicator can undergo a color change in a substantial temperature range, namely from a very low temperature (e.g., -30°C) to a very high temperature of a few hundred degrees centigrade. The preferred temperature is below about 60°C, most preferred temperature is ambient.

<u>Time</u>: The time required for color change depends upon the nature of the plasma. A more reactive and more diffusive plasma would change color faster. The time required for the color change can be varied by controlling the diffusion of the components of the plasma by adding proper additives in the polymer. Crosslinking the polymer will reduce the diffusion of the plasma and increase the time required for the color change. Additives which react with plasma can also increase the time required for the color change. Similarly, certain additives, such as oxidant and plasticizer would decrease the time

25

19

required for the color change. The time required for the color change can also be varied by selecting proper indicators, activators, polymers and their mixtures.

For a given sterilization cycle, the time required for the color change

10 can be varied by varying one or more of the following parameters:

w/w = weight to weight percent

1. Thickness of the polymer indicator layer.

The thickness of the indicator and barrier layers may vary from a micron to a few hundred microns. The preferred thickness is 1-50 microns and the most preferred range is 2-20 microns.

2. Concentration of the activator.

The concentration of activator may vary from 0.1 to 50 w/w%. The preferred concentration is 1 to 20 w/w% and the most preferred concentration is 2-10 w/w/%.

3. Concentration of the indicator.

The concentration of the indicator may vary from 0.1 to 20w/w%. The preferred concentration is 10 to 10w/w% and the most preferred concentration is 2-5 w/w%.

4. Concentration of other additives.

The concentration of additives may vary from 0.1 to 20w/w%. The preferred concentration is 0.5 to 10w/w% and the most preferred concentration is 1-5w/w/%.

5. Nature of the polymer

35 6. Nature of the barrier

This is the same group as the polymer, but may be different from it

7. Thickness of the barrier

The barrier may be between 2 and 200 microns thick, preferrably 2-20 microns.

10 8. Nature of the activator

9. Nature of the indicator

10. Nature of the additives

15

he preferred time range for sterilization is from 5 minutes to a few hours, The most preferred time is about 15-60 minutes.

In addition to the plasma of hydrogen peroxide, the device may be used with other plasmas, such as that of peracetic acid and a mixture of hydrogen, oxygen and argon, as well as for hydrogen peroxide itself.

<u>Advantages</u>: The plasma indicator disclosed here offers the following advantages:

- 25 *. It is selective to plasma (i.e., no or least effect of other sterilants).
 - *. It provides desired color changes (from a starting light color, such as white/colorless, yellow, orange, pink, or red to a final dark color, such as blue, green, black, purple or violet).
 - *. It provides an intermediate color for monitoring a partial cycle.
- *. The time, temperature etc required for the color change can be varied by simple means.
 - *. There is essentially no effect of ambient conditions (e.g., dry heat, humidity and light) before and after the sterilization.
 - *. It is unaffected by sealing hot bar.
- 35 *. It has required pot life.
 - *. There is no bleeding/diffusion of dyes.

- *. The ingredients (indicators/dyes and activators/additives) are water soluble.No grinding of ingredients required.
 - *. Formulations can be made by simple procedures (mixing/dissolution).
 - *. Ink is printable with gravure and flexo presses on polyester, paper and Tyvek.
- 10 *. The print rolls are easy to clean.
 - *. It uses least toxic or hazardous chemicals.
 - *. It uses readily available chemicals (dyes, activators and binders).
 - *. There is an option of precision measurement with moving boundary.
 - *. Formulations are inexpensive.
- 15 *. It is unaffected by ethylene oxide, steam, heat and radiation.

5

EXAMPLES

Example 1. General procedure for preparation of the sample devices.

10

15

In a 10ml test tube were added 3 ml of EC001270 (an acrylate ink extender supplied by Environmental Inks and Coating, Co, Lithicum, MD), about 0.5 ml of solution (60 w/w% TEAB in water) and 0.5 ml of a solution (e.g., 4 w/w% solution of a phenol red in ethanol). The contents were mixed and coated with a #10 wire wound rod on a 100 micron polyester film, paper and Tyvek^R. The coatings were dried in an oven at about 50°C for 5 minutes.

Example 2. Exposure to vapor of hydrogen peroxide.

20

Samples of example 1 were hung on the side of a 20-liter clear glass battery jar containing an empty petri dish. About 50 ml of 30% hydrogen peroxide was poured in the petri dish and the jar was tightly closed. The color changes of the samples were noted. For testing under low pressure a vacuum desicator was used instead of the battery jar.

25

Example 3: Effect of activator on dyes.

Using the general procedure described in example 1, coatings were prepared from EC001270 as a polymer, TEAB as an activator, and several of the dyes listed in Table 2 as indicators. The coatings were exposed to one sterilization cycle of hydrogen peroxide plasma using STERRAD® model 100, sterilization system made by Advanced Sterilization Products (address). Most of these dyes changed their colors. Some of them either became lighter or colorless while the others developed colors. Only a small number of the dyes showed color change without the activator. Representative color changes are listed Table 4.

35

5 .Table 4. Representative color changes of some dyes with ECOO1270 and TEAB

	Dye	Original color	Color after plasma
			treatment
10	Acid blue 93	Blue	Light blue
	Evans blue	Blue	Coloriess
	Thionin	Blue	Red
	Disperse blue 3	Blue	Light red
	Malachite green oxa	alate Green	Coloriess
15	Disperse red 14	Purple	Yellow
	Acid fuschin	Red	Coloriess
	Pararosaniline acet	ate Pink	Colorless
	Dichlorofluorescein	Light pink	Red
	Brilliant yellow	Yellow	Red
20	m-Cresol purple	Yellow	Purple
	Phenol red	Yel low	Blue
	Cresol red	Yellow	Green
	Fluorexon	Yellow	Red
	Basic yellow 11	Yellow	Colorless
25	Fluorescein	Coloriess	Red

Similar color changes were obtained with many other activators, such as sodium nitrite, sodium thiocyanate, salts of amines and acids and other binders such as polyurethane (Witcobond 213, Witco Corporation, Houston,

30 TX) and cellulose nitrate.

Example 4. Effect of other halides.

Using the procedure described in example 1, coatings were prepared with certain organic and inorganic halides as activators and fluorescein, phenol red, chlorophenol red, acid fuschin, brilliant green and their mixtures as indicators. The coatings were exposed to hydrogen peroxide and its plasma. All coatings of acid fuschin changed from red-to-colorless and that of brilliant green changed from green-to-colorless with plasma. The color changes of the other three dyes are listed in Table 5.

<u>Table 5</u>

<u>Effect of different halides on color changes of some dyes with hydrogen</u>

<u>peroxide plasma.</u>

20 The original color of the coatings was yellow

Dye	Fluorescein	Phenol Red	Chlorophenol Red
Acetylcholine chloride	Pink	L.blue	L.Yellow-green
Ammonium bromide	Pink	Green	L.Yellow-green
Choline chloride	Orange	Colorless	L.Yellow-green
Tetrabutyl ammonium bromide	Pink	Yellow	L.Yellow-green
Choline iodide	Orange	Colorless	L.green
Potassium bromide	Pink	Blue	L.green
Potassium iodide	Yellow	Yellow	Yellow
Sodium iodide	Yellow	Yellow	Yellow
Tetraethylammonium bromide	Pink	Blue	Blue green
Tetrahexylammonium bromide	Yellow	Yellow	Yellow
Tetramethyl ammonium acetate	e Orange	Blue	L.blue

L. = light

Mixtures of halides also provided essentially the same color changes. Bromides were more effective than other halides in introducing the color

change. The iodides show color change when polyvinyl alcohol is used as a 5 binder.

Example 5. Effect of bromides.

Using the general procedure described in example 1, coatings were prepared using phenol red as an indicator and several organic and inorganic bromides e.g., ammonium bromide, barium bromide, calcium bromide, iron (III) bromide, potassium bromide, zinc bromide, diphenyliodinium bromide, tetraethylammonium bromide, bromide. tetrabutylammonium tetrahexylammonium bromide, dimethyldioctadecylammonium bromide and 15 phenacylpyridinium bromide and mixtures of some bromides as activators. The coatings were exposed to hydrogen peroxide vapor and its plasma.

All coatings changed from yellow-to-blue. Coatings containing dimethyldioctadecylammonium bromide changed slowly while those of ammonium bromide, potassium bromide and tetraethylammonium bromide changed faster.

Example 6. Effect of pH of the polymer medium.

25

30

20

10

The pH of the EC001270 coating is about 4-5 while that of polyvinyl alcohol (PVOH) is about 7. Coating were prepared using potassium iodide and TEAB as activators and brilliant yellow, m-cresol purple, cresol red and phenol red as indicators in EC001270 and polyvinylalcohol as binders. The The coatings of coatings were exposed to hydrogen peroxide vapor. KI/PVOH changed colors in seconds. The results indicate that the device can be used for monitoring hydrogen peroxide. As shown in Table 6, depending upon the pH of the binders, different color can be obtained.

Table 6. Color of some dyes with different polymer and activators.

10	Dye	Polymer Activator	EC00127	70	PVOH TEAB	EC001270 KI	PVOH KI	
	Brilliant y m-Creso Cresol R	ol purple	Light red Blue Light blu		Yellow Light orange Yellow	Light orange Orange Blue	Red Blue Blue	
15	Phenol F	Red	Blue	Yellow	Yel	low	Purple	

Effect of nature of polymer: The other polymers were much slower, e.g., took more than ten hours compared to less than hour for EC.

20

30

Example 7. Effect of a barrier coat.

A coating of EC001270, phenol red and TEAB bromide of example 1 was top-coated with EC001270 using number 3, 5, and 10 wire wound rods. After drying, the coatings were exposed to hydrogen peroxide vapor. The time required for the color change is shown in the table 7.

Table 7. Effect of thickness of a barrier coat on the time required for the color change.

	Top-coating Bar	Time (hours) required for the color change
35	None	0.1
	3	3
	5	7
	10	18

10

Example 8. Two layer device, both containing indicators.

coating of EC001270, phenol red and TEAB was top-coated with that of EC001270, pararosaniline acetate and TEAB. This double-coated device appeared red. The device was exposed to hydrogen peroxide vapor. Pararosaniline of the top-coat changed from red-to-colorless. Phenol red of the undercoat then changed from yellow-to-blue. The color changes of the device with time are shown in Table 8.

15 <u>Table 8. Color changes of the two-layer device having phenol red in the undercoat and pararosaniline in the top-coat.</u>

Time (hours)	0	3	6	10	24
Color	Red	Light red	Purple	Blue	Dark blue

Example 9. Two-layer device, one containing indicator.

25

20

The two layer device was prepared by first coating a red pigmented ink followed by a top-coating of EC001270, malachite green oxalate and TEAB. The resultant device appeared very dark, almost black. The device was exposed to hydrogen peroxide. The device changed from black-to-red in about six hours. The color change is due to fading of the green color of malachite green oxalate which make the red color of the under coat visible.

Example 10. Moving boundary device.

35

30

A moving boundary device was created by applying a wedge shaped coating with a wedge shaped coating bar (0 at one end and #10 at the other

end) of EC001270 on the coating of EC001270, phenol red and TEAB. The device thus created was exposed to hydrogen peroxide vapor. A boundary was created between the original yellow color and the new blue color after about 2 hours. The boundary moved from the thin end towards the thick end. The distance traveled by the boundary is shown in Table 9.

10

Table 9. Distance traveled by the moving boundary device

	======================================						
	Time (hours)	2	4	10	24		
15	Distance (cm)	0.2	0.5	2	3		

Example 11. Detection of residual absorbed hydrogen peroxide.

A piece of 100-micron polyester film having a coating of EC001270 was exposed to hydrogen peroxide vapor for two hours. The coated film was removed and hung in air for ten minutes. A coating of polyvinyl alcohol, m-cresol red, and KI coated on polyester film was placed on hydrogen peroxide exposed film. The yellow coating turned light blue within 15 minutes.

In accordance with the above procedure, the device may be used for the detection of ozone and other oxidative vapors.

5 Example 12. Pilot coating.

In 1200g EC001270 extender the following were added while stirring: 18g of 28% ammonium hydroxide, 100g of water, 200ml of isopropanol, 36g pararosaniline acetate, 18g of phenol red and 64g of TEAB. The formulations were coated on Tyvek^R and polyester film using a pilot coater of Rexam Medical Packaging, Mt. Holly, NJ. Eleven other formulations with different ratios of the dyes, TEAB and crosslinking agent (a polyaziridine) were also coated. The coatings were red (burgundy) color.

The samples were treated with half, one and two cycles of plasma sterilization cycles of hydrogen peroxide plasma using STERRAD® model 100, sterilization system of Advanced Sterilization Products. The samples treated with half cycle were yellow-green while those treated with one and two cycles were blue.

20

25

30

35

10

15

The coated samples were exposed to 100% ethylene oxide and steam sterilization cycle. We found that there was no or negligible effect of (1) 100% ETO for 16 hours and (2) normal steam sterilization treatment (e.g., 20 minutes at 121°C). Steam makes the color a bit lighter and ETO makes slightly purple/violet.

The samples were annealed at 60°C for five days and exposed to 100% humidity at 60°C for 5 days and then treated with the plasma along with control (un-treated samples) for determination of the shelf life and color stability. We found that there was no or negligible effect of: (1) dry heat at 60°C five days and (2) 100% humidity at 60°C for five days.

Samples treated with half, one and two cycles of the plasma were annealed at 60°C for five days and exposed to 100% humidity at 60°C for 5 days for determination archival life and color stability. There was no or negligible effect of (1) dry heat at 60°C five days, (2) 100% humidity at 60°C for five days on the final blue color.

The effect of hydrogen peroxide vapor depends on the exposure time. There was essentially no color change for about 4 hours. After about 5 hours the samples started turning yellow. At about 8 hours it started turning green color and at about 16 hours, they turned blue.

10

Example 13: Alternate embodiment

A coating of EC001270, phenol red and TEAB was top-coated with that of EC001270, pararosaniline acetate and TEAB. This double-coated device appeared red. The device was exposed to hydrogen peroxide vapor and its plasma. First, pararosaniline in the top-coat changed from red-to-colorless. Phenol red in the undercoat then changed from yellow-to-blue.

20

25

30

15

Example 14: Alternate embodiment

An alternate embodiment of this type of device was prepared by first coating a red pigmented ink followed by a top-coating of EC001270, malachite green oxalate and TEAB. The resultant device appeared very dark, almost black. The device changed from black-to-red when exposed to hydrogen peroxide.

To get orange-to-yellow and red/burgundy-to-yellow, methyl red was mixed with cresol red and acid fuschin with cresol red along with TEAB as activator and EC001270 as a polymer. Red-to-blue color was obtained by mixing pararoseaniline acetate and phenol red as shown in example 12.

Example 15: Alternate embodiment

A mixture of pararoseaniline acetate and phenol red in the EC extender with tetraethylammonium bromide as an activator provides

5 burgundy-to-green-to-blue color changes when treated with hydrogen peroxide plasma.

To get orange-to-yellow and red/burgundy-to-yellow, methyl red was mixed with cresol red and acid fuschin with cresol red along with TEAB as activator and EC001270 as a polymer.

While this invention has been shown and described in detail in the context of a preferred embodiment, and with various modifications thereto, a wide variety of other modifications can be made without departing from scope of the inventive teachings.

AMENDED CLAIMS

[received by the International Bureau on 2 October 2000 (02.10.00); original claims 1, 40 and 41 amended; remaining claims unchanged (2 pages)]

5

CLAIMS

I claim:

- 10 1. A device for monitoring an oxidizing vapor or plasma comprising: at least one layer of polymer, having incorporated therein
 - a) an indicator capable of undergoing at least one color change
 - b) an activator for said indicator wherein said activator, when contacted with said oxidizing plasma, undergoes a reaction wherein the product of said reaction causes said indicator to undergo said color change.
 - 2. The device of claim 1 wherein the said indicator comprises at least one member of the group consisting of pigments, dyes, precursors of said dyes, and mixtures of any of the foregoing group members.

- The device of claim 1 wherein the said indicator is a pH-sensitive sensitive dye.
- 4. The device of claim 1 wherein the said indicator is phenol red, π 25 cresol purple, pararosaniline or mixtures thereof.
 - 5. The device of claim 1 where the said indicator undergoes halogenation or oxidation.
- 30 6. The device of claim 1 wherein the said indicator undergoes a yellow-to-blue, red-to-yellow or red-to-blue color change.
 - 7. The device of claim 1 wherein said polymer is soluble in an organic solvent.

15

25

- 5 8. The device of claim 1 wherein said polymer is soluble in water or is water dispersible.
 - 9. The device of claim 8 wherein said polymer is a water soluble or water dispersable homopolymer, or a copolymer or a mixture thereof.
 - 10. The device of Claim 1 wherein said polymer is a polymer of styrene, acrylate, acrylic acid, acrylamide, vinyl acetate, vinyl alcohol, vinyl chloride, styrene, polyurethanes, cellulose nitrate, carboxymethyl cellulose or a mixture thereof.
 - 11. The device of claim 10 wherein said polymer is a homopolymer, or a copolymers or a mixture thereof.
- 12. The device of claim 8 wherein said polymer is a polymer of styrene, acrylate, acrylic acid, acrylamide, vinyl acetate, vinyl alcohol, vinyl chloride, styrene, polyurethanes, cellulose nitrate, carboxymethyl cellulose or a mixture thereof.
 - 13. The device of claim 1 wherein the polymer is an acrylate polymer.
 - 14. The device of claim 1 wherein the polymer is cellulose nitrate or carboxymethylcellulose
- 15. The device of claim 1 wherein the reaction product of said activator30 and said plasma is a halo-acid.
 - 16. The device of claim 1 wherein the said activator is a salt.
 - 17. The device of Claim 1 wherein said activator a halide.
 - 18. The device of Claim 1 wherein said activator is a bromide.

- 19. The device of claim 1 wherein the said activator is a bromide of alkali metal or quaternary amine.
- 20. The device of claim 1 wherein said activator is tetrabutylammonium bromide or tetraethylammonium bromide or mixture thereof.
 - 21. The device of claim 1 wherein said activator is a salt of an amine and an organic or inorganic acid acid
- 15 22. The device of claim 1 wherein said activator is a thiocyanate
 - 23. The device of claim 1 wherein said activator is sodium thiocyanate
- 24. The device of claim 1 additionally comprising an additive to control the diffusion of plasma gases.
 - 25. The device of claim 1 additionally comprising a crosslinking agent or a plasticizer to control the diffusion of plasma gases.
- 25 26. The device of claim 1 additionally comprising a zinc compound or a polyaziridine to control the diffusion of plasma gases.
 - 27. The device of claim 1 comprising two layers.
- 30 28. The device of claim 1 additionally comprising a polymeric top layer.
 - 29. The device of claim 1 additionally comprising a wedge shaped polymeric top layer.
- 35 30. The process of making a device of claim 1 which comprises dissolving or dispersing the components thereof in a solvent therefor, applying the thus

- formed solution or dispersate to a substrate and permitting the solvent to evaporate.
 - 31. The process of claim 30 wherein the substrate is a container for an item to be sterilized.

- 32. The process of claim 30 wherein the substrate is a plastic film, paper or metal.
- 33. The process of claim 30 wherein the substrate is polyester film or spunbonded polyolefins.
 - 34. The process of claim 30 wherein the solution is an ink formulation.
- 35. The process of claim 30 wherein the solution is an aqueous ink 20 formulation.
 - 36. The process of claim 35 said ink formulation comprises an acrylate polymer.

25

- 37. A process of using a device of claim 1 for monitoring sterilization of materials comprising the steps of
- a) affixing the device to said materials or containers containing same
- b) carrying out the process of sterilization including the step of introducing
 30 the plasma into a vessel containing said materials or containers therefore and
 - c) observing the presence of a color change of said device.
 - 38. The process of claim 37 wherein the plasma is derived from a member selected from the group consisting of hydrogen peroxide, perchloric acid and oxygen.

- 5 39. The process of claim 37 wherein the plasma is that derived from hydrogen peroxide
 - 40. A process of using the device of claim 1 for monitoring an oxidizing vapor comprising the steps of
 - a) exposing the device to an oxidizing plasma,
 - b) observing the presence of color change in the device.
 - 41. The process of claim 40, wherein the oxidizing plasma is ozone or hydrogen peroxide.




Figure 1




Figure 2




Figure 3






Figure 4

DECLARATION F	אר וודוו פא	Alterney Docket Num	per PA	rl3.0-010		
DESI	ON	First Named Inventor				
PATENT AP	3, 1	COMPLETE IF KNOWN				
(37 CFR	(1.63)	Application Number				
Declaration	Declaration	Filing Date				
Submitted OR	Submitted after Initial Filling (surcharge	Art Unit				
with Initial Filing	(37 CFR 1.16 (e)); required)	Examiner Name				
the below named inventor, I						
rveiden æ, melling address, an						
leve I I'm the original and first	inversor of the subject matter	which is claimed and for which	8 palent is so	oughl on the invention and		
	والمراجع والمجاهد المراس الماسك والمستكاد المسائد	and the many of the state of th				
dicators for r	monitoring ste	rilization wit	n Plasm	ıa		
	والمحق أتعلم الم	a single				
	er en	The state of the s				
	*					
		The state of the s				
	774- A-1-	6.5.5				
	(Title of the	(אסממשייוו				
specification of which		``,				
is akached hersto	0					
IN GRAPHAN LINERED	4.1	Siere Charles				
OR		18 2 21				
was fled on (MM/DD/YYYY)	04 (33 (3000	pt Johad States and	lastica Al	nom SVT International		
not an family () ()	04/11/2000	. BL United States App		www.relinumauなる		
			(0.0.40			
Callon Number PCT/US	007094 grd was amend	an (MM/DD/YYY) 10/	02/200	(if applicable).		
Carried Ma						
by state that I have reviewed a	and understand the contents of	of the above identified anenifica	tion including	the claims as amended i		
mendment specifically referred	to above.			in the mention of the Child (1990)		
towledge the duty to disclose is	Mortration which is material h	naterialist pe defined in 27	CER 1 56 inc	duding for continuetion is		
iowleage the duty to disclose it ations, material information wh ational fiting date of the contin	ich became evallable betwee	in the filing date of the prior ap	cation and t	he national or PCT		
MANAGEMENT CHEST OF CHEST CONTINUE	auon-in-part application					
by deim foreign priority benefit ers rights certificate(s), or 36	is under 35 U.S.C. 119(a)-(d)	or (I), or 365(b) of any toroign	application(s) for palent, inventor's or p		
COLANTIACION HATALINAMENTALINA	DAVE SIGN MEMBER PRINT.	by chartico the hov any formi	- acclication	for national lationality or a		
er's nghis certificate(s), or an	y PCT International applicati	on having a filing date before	that of the a	polication on which prior		
hiar Foreign Application	7	Korelno Sili D-4	Dalante.	Contrad Command		
Number(s)	Country	Foreign Filing Date (MM/DD/YYYY)	Priority lot Claimed	Certified Copy Attach		
· · · · · · · · · · · · · · · · · · ·	4			YES NO		
		,	ليا			
		1	ليسييا			
	,	· ·				
	- 101					
		1				

DECLARATION — Utility or Design Patent Application

Direct all correspondence to: Customer Number or Bar Code Label OR Correspondence address below						
	01815					
Name	PATENT TRADEMARK OF	FICE				
	at Washington					
Address						
Address						
City	State		ZIP			
City	State		ZIF			
Country	elephone		Fax			
I hereby declare that all statements made herein of my are believed to be true; and further that these stateme made are punishable by fine or imprisonment, or both, validity of the application or any patent issued thereon.	nts were made with the k	nowledge that willful false	statements and the like so			
NAME OF SOLE OR FIRST INVENTOR :	A petition has be	en filed for this unsigr	ned inventor			
Given Name (first and middle [if any]) Gordhanbhai N or Surname Patel						
Inventor's Signature Date						
Residence: City Middlesex	State NJ	USA Country	US Citizenship			
Mailing Address 120, Wood Avenue	3					
_{City} Middlesex	State NJ	2IP 08846	Country USA			
NAME OF SECOND INVENTOR:	A petition has been	n filed for this unsigne	d inventor			
Given Name (first and middle [if any]) Family Name or Surname						
lavoré de la companya						
Inventor's Signature			Date			
Residence: City	State	Country	Citizenship			
	. , , , , , , , , , , , , , , , , , , ,	•	'			
Mailing Address						
City	State	ZIP	Country			
	upplemental Additional Inv	entor(s) sheet(s) PTO/SB/	·			

PTO/SA/01 (10-01)

Approved for use through 10/3 U20/2. Dick 8 8841-0032

U.S. Parent and Trademark Celice; U.S. DEPARTMENT OF COMMERCE

U.S. Parent and Trademark Celice; U.S. DEPARTMENT OF COMMERCE

Under the Paperwork Reduction Act of 1895, no persons are required to respond to a collection of information unless it contains a ward CMS control number

DECLARATION — Utility or Design Patent Application Direct eli correspondence to: Customer Number or Bar Code Lubel OR Correspondence address below 01815 Name PATENT TRADSMARK OFFICE City ZIP Country Telephone I hereby designe that all statements made herein of my own knowledge are true and that all statements made on information and belief are belief for the statements made with the knowledge that wilful false statements and the second made are pt. schable by fine or impresement, or both, under 18 U.S.C. 1001 and that such willful false statements may jeopart to the validity of the application or any patent issued thereon. NAME OF SOLE OR FIRST INVENTOR: A patition has been filed for this unsigned inventor Family Name (first and middle (Many)) Gordhanbhai or Surname **Patel** Data 14/01 Ford NunWhee Middlesex NJ US Residence: Fity. 120, Wood Avenue Mailing Adda one ZP 08846 Middlesex , NJ USA NAME OF SECOND INVENTOR: *** A petition has been filed for this unsigned inventor GMen Name Family Name (first and micide (if any)) or Sumame inventor's Signature Residence: (Sily Ckizenehlp Mailing Addrass Chy Country supplemental Additional Inventor(s) sheet(s) PTO/SB/12A attected hereto. Additional inventors are being ramed on the

[Page 2 of 2]

now how in

United States Patent & Trademark Office

Office of Initial Patent Examination -- Scanning Division



Application deficiencies found during scanning:

Page(s) of were not present for scanning.

Document title)

Page(s) of were not present for scanning.

Scanned copy is best available. Ligures are dark.